



### Sample requirements for DNA sequence capture libraries

NimbleGen SeqCap EZ DNA sequence capture libraries.....	Page 2
Agilent SureSelect XT DNA sequence capture libraries.....	Page 3
Agilent SureSelect XT HS and Low Input DNA sequence capture libraries.....	Page 4




### Sample requirements for NimbleGen SeqCap EZ DNA sequence capture libraries

- The following amounts of **unamplified** genomic DNA are required for each sequence capture:

Sequencing platform	Capture method	DNA quantity
Illumina	NimbleGen SeqCap EZ - Standard input	2000 ng in 50-100 µl of 1xTE buffer
Illumina	NimbleGen SeqCap EZ - Low input	200 ng in 50-100 µl of 1xTE buffer

- DNA concentration should be determined by a fluorometry-based assay such as Qubit or PicoGreen.
- The  $A_{260}/A_{280}$  and  $A_{260}/A_{230}$  ratios (which are obtained by NanoDrop) must be  $\geq 1.8$  and  $\geq 1.9$ , respectively, for all samples (although this may not be possible for the samples that are to be processed through the low input protocol).
- The DNA sample should not contain particulate matter.
- The DNA must be high quality, high molecular weight and free from RNA contamination. This should be assessed by electrophoresis of 50-100 ng of each sample on a 1% agarose gel alongside molecular markers of the appropriate size. Genomic DNA should be visible as a single prominent band  $>12$  Kb. If there is more than one band or a smear, the DNA may be degraded or contain a contaminant that could affect the library preparation. RNA will appear as a smear  $<200$  bp. If RNA is present, please treat the DNA with RNase A (DNase-free) and then purify, preferably using a column based method.

The gel image below illustrates intact, degraded and RNA contaminated genomic DNA run alongside Hyperladder I (Bioline). If a smear is in the range between zero and  $\sim 2000$  bp, RNase treatment may be required to reveal whether the degraded material is gDNA or RNA.



**For projects that involve  $\geq 32$  samples**, please submit samples in a well sealed 96-well plate to match our workflow.

**For projects that involve  $\leq 31$  samples**, please label the samples and tubes 1, 2, 3 etc. in order to ease our sample identification. Please remember to underscore numbers that can be read upside-down such as 6, 9, 16, 91, 69, 96 etc.




### Sample requirements for Agilent SureSelect XT DNA sequence capture libraries

- The following amounts of **unamplified** genomic DNA are required for each sequence capture:

Sequencing platform	Capture method	DNA quantity
Illumina	Agilent SureSelect XT- Standard input	3000 ng in 50-100 µl of 1xTE buffer
Illumina	Agilent SureSelect XT - Low input	200 ng in 50-100 µl of 1xTE buffer

- DNA concentration should be determined by a fluorometry-based assay such as Qubit or PicoGreen.
- The  $A_{260}/A_{280}$  and  $A_{260}/A_{230}$  ratios (which are obtained by NanoDrop) must be  $\geq 1.8$  and  $\geq 1.9$ , respectively, for all samples (although this may not be possible for the samples that are to be processed through the low input protocol).
- The DNA sample should not contain particulate matter.
- The DNA must be high quality, high molecular weight and free from RNA contamination. This should be assessed by electrophoresis of 50-100 ng of each sample on a 1% agarose gel alongside molecular markers of the appropriate size. Genomic DNA should be visible as a single prominent band  $>12$  Kb. If there is more than one band or a smear, the DNA may be degraded or contain a contaminant that could affect the library preparation. RNA will appear as a smear  $<200$  bp. If RNA is present, please treat the DNA with RNase A (DNase-free) and then purify, preferably using a column based method.

The gel image below illustrates intact, degraded and RNA contaminated genomic DNA run alongside Hyperladder I (Bioline). If a smear is in the range between zero and  $\sim 2000$  bp, RNase treatment may be required to reveal whether the degraded material is gDNA or RNA.



**For projects that involve  $\geq 32$  samples**, please submit samples in a well sealed 96-well plate to match our workflow.

**For projects that involve  $\leq 31$  samples**, please label the samples and tubes 1, 2, 3 etc. in order to ease our sample identification. Please remember to underscore numbers that can be read upside-down such as 6, 9, 16, 91, 69, 96 etc.



### Sample requirements for Agilent SureSelect XT HS and Low Input DNA sequence capture libraries

- The following amounts of **unamplified** genomic DNA are required for each sequence capture:

Sequencing platform	DNA sample source	DNA quantity
Illumina	Fresh biological samples (non-FFPE)	10-200 ng in 30-60 µl of buffer AE*
Illumina	FFPE tissue sections	All that is available in 30-60 µl of buffer ATE**

\* Buffer AE composition:

10 mM Tris-Cl

0.5 mM EDTA; pH 9

\*\* Buffer ATE composition:

10 mM Tris-Cl, pH 8.3

0.1 mM EDTA

0.04% NaN<sub>3</sub>

- For fresh biological samples (non-FFPE), please refer to page 3 for sample QC requirements.
- For FFPE samples, DNA concentration should be determined by a fluorometry-based assay such as Qubit or PicoGreen. A qPCR-based assay will be used to quantify and qualify these samples at the CGR prior to library preparation.