

---

3T MRI Scanning Procedure  
Study: CNS-COVID

---

**Author:**      **Name**                      Jon Campbell  
                         **Position**                      Lead Research Radiographer  
                                                              Wellcome Centre for Integrative Neuroimaging  
                                                              University of Oxford

**Checked:**    **Name**                                  Peter Jezzard  
                         **Position**                                  Professor of Neuroimaging  
                            Wellcome Centre for Integrative Neuroimaging  
                            University of Oxford

### Document History

Version	Initials	Date	Comment
1	JC	04/01/21	

### Related Documents and Location

PIS and Consent Form (electronic, study docs folder)

---

## Study Information

**Study Title:** multi-modal magnetic resonance imaging (MRI) protocols at 3 Tesla for use in multi-site COVID-19 research

**Ethics Number:** Cambridge: HBREC.2020.44

**Scanner:** 3T

**Session Time:** 60min

---

## Equipment Required

### Standard

- 32 Channel Head Coil
- Regular mirror
- Ear plugs
- Immobilization sponges
- Buzzer

### Lighting and stimulus

- Bore lights on
- Room lights on
- LCD on (Biobank fixation cross)

### Ancillary

- Siemens Respiratory bellows (connects to the ECG unit)
- Siemens Pulse meter

[These should be used if available, but are not mandatory – if used the physiological information during the DTI and fMRI scans is saved to the DICOM files]

## Protocol Sequences

### 1. -- Centre to the coil centre-mark --

Note that all sequences are set to REF or FIX with a 0mm offset, such that the table does not move between scans (to avoid re-shimming).

### 2. AAHead\_Scout


14s

- Auto-runs, no setup required
- Once finished sagittal, coronal and axial reconstructions will automatically load into the graphic segments (as below)



### 3. T1\_p2\_1mm\_fov256\_sag\_TI\_880

4m54s



- Open scan, allow auto-align to position scan, click close and run scan
- NB scan will automatically shim


### 4. -- fixation cross --

[Display the Biobank fixation cross during the resting state fMRI scan and ask the subject to try not to think about anything specific, to stay awake with eyes open, and to fixate on the cross. The fixation cross can be downloaded from:

<https://www.fmrib.ox.ac.uk/ukbiobank/protocol/Crosshair.png>]


### 5. -- Shim x3 --

[Perform a manual scan three times to obtain a good starting shim over the brain. To do this, open the MB8\_FMRI\_fov210\_2.4mm\_resting sequence and click on 'Options -> Adjustments' to bring up the manual shimming interface



- Now click on 'Options' -> 'Adjustments' to bring up the manual shimming interface

- Select the 3D Shim tab (yellow arrow below)



- Shim by clicking on 1) Measure, 2) Apply, 3) Calculate and 4) Apply (as indicated above).

- Repeat this 2 more times and you should see 6 rows of text in the light grey shim box, along the lines of:


## Manual Adjustments

No	X [ $\mu\text{T/m}$ ]	Y [ $\mu\text{T/m}$ ]	Z [ $\mu\text{T/m}$ ]	A20 [ $\mu\text{T/m}^2$ ]	A21 [ $\mu\text{T/m}^2$ ]	B21 [ $\mu\text{T/m}^2$ ]	A22 [ $\mu\text{T/m}^2$ ]	B22 [ $\mu\text{T/m}^2$ ]	Converged
1	-	-	-	-	-	-	-	-	Yes
2	792.69	-1,847.40	-216.32	286.25	41.43	6.92	-36.98	4.08	Yes
3	-	-	-	-	-	-	-	-	Yes
4	792.32	-1,847.79	-215.85	292.68	40.52	11.45	-37.98	1.48	Yes
5	-	-	-	-	-	-	-	-	Yes
6	792.56	-1,847.52	-215.88	291.16	42.38	13.65	-38.44	3.74	Yes

- Close and apply to start the scan
- Remind researcher to display the fixation cross
- Resting state instructions: **“Look at the cross for the next scan, blink normally and try not to fall asleep”**


## 6. MB8\_FMRI\_fov210\_2.4mm\_resting

6m10s



## 7. diff\_PA\_MPopt\_MB3\_3b0\_lowflip


36s



- Open scan, allow auto-align to position scan, click close and run scan
- NB do not change phase encoding direction

### 8. diff\_AP\_MPopt\_MB3\_50b1000\_50b2000\_8b0\_lowflip


6m32s



- Open scan, allow auto-align to position scan, click close and run scan
- NB do not change phase encoding direction

### 9. t2\_space\_dark-fluid\_sag\_p3

4m32s



- Open scan, allow auto-align to position scan, click close and run scan

### 10. SWI\_3mm\_Updated\_v1.1


2m8s



- Open scan, allow auto-align to position scan, click close and run scan

### 11. fme\_pCASL\_M0\_RL


16s



- Open scan, allow auto-align to position scan, click close and run scan

### 12. fme\_pCASL\_BL1800\_PLD400

26s



- Open scan, allow auto-align to position scan, click close and run scan

### 13. fme\_pCASL\_BL1800\_PLD800

30s

- Open scan, allow auto-align to position scan, click close and run scan
- NB appearance identical to **fme\_pCASL\_BL1800\_PLD400**

### 14. fme\_pCASL\_BL1800\_PLD1200

34s

- Open scan, allow auto-align to position scan, click close and run scan
- NB appearance identical to **fme\_pCASL\_BL1800\_PLD400**

### 15. fme\_pCASL\_BL1800\_PLD1600

38s

- Open scan, allow auto-align to position scan, click close and run scan
- NB appearance identical to **fme\_pCASL\_BL1800\_PLD400**

### 16. fme\_pCASL\_BL1800\_PLD2000

42s

- Open scan, allow auto-align to position scan, click close and run scan
- NB appearance identical to **fme\_pCASL\_BL1800\_PLD400**

### 17. -- Do not change angle --

- Reminder to not change the default angle or straighten the following 2 sequences

### 18. -- Manually position so labelling


### 19. plane 2cm below cerebellum --

### 20. -- Click OK on conflict prompt --

- The `fme_pCASL_GE_M0` has a known conflict prompt which can be ok'd

## 21. `fme_pCASL_GE_PLD2025`

4m46s




- This scan must be manually positioned
- ! Do not change angle or straighten this sequence !
- Move slice slab so that labelling plane is 2cm below the cerebellum and click apply
- NB when you click 'Apply' a yellow triangle will appear on the following sequence. This can be ignored

17	<code>fme_pCASL_GE_PLD2025</code>	17	04:46
18	<code>fme_pCASL_GE_M0</code>	17	00:50

## 22. `fme_pCASL_GE_M0`

50s



- Auto-copies COSG&SR from sequence `fme_pCASL_GE_PLD2025`
- Click OK on the conflict prompt when it appears

