

Looking after your prism

(Easy read version)

Cleaning your prism

The rough side can get dusty and will need cleaning.

- Use cool water slowly from a tap and a soft toothbrush – brush with the lines on the prism.
- Pat it dry with a soft cloth
- If very dirty, remove and clean with a drop of washing-up liquid and water.
- Rinse with water



Re-fitting your prism

If you need to put your prism back on

- Put your glasses and the prism under a slow running tap of cool water (1).
- Slide smooth side onto the back of your glasses lens (2).
- Take your glasses out of the water, holding the prism in place.
- Gently squeeze out any bubbles with your thumb (3).
- The prism may slide around now if you need to change where it is.
- Pat dry using a clean, dry paper towel.
- Be careful with your glasses until the prism has dried on.

